

CSI ClorDiSys Solutions, Inc.™

"The Chlorine Dioxide People"

What is Chlorine Dioxide? Where is it used? How does it work?

General

Chlorine Dioxide (CD) is a greenish-yellow gas and is a single-electron-transfer oxidizing agent with a chlorine-like odor. CD has been recognized since the beginning of the century for its disinfecting properties; and has been approved by the US EPA for many applications including the widespread use of CD in the treatment of drinking water. Beyond this and numerous other aqueous applications, the sporicidal properties of *gaseous* CD were demonstrated in 1986. Subsequent to these initial studies, it has been shown that gaseous CD is a rapid and effective sterilant active against bacteria, yeasts, molds, and viruses. The rapid sterilizing activity of CD is present at ambient temperature and at relatively low gas concentration, 1 to 30 mg/L.

Uses

Chlorine dioxide is widely used as an antimicrobial and as an oxidizing agent in drinking water; poultry process water, swimming pools, and mouthwash preparations. It is used to sanitize fruit and vegetables as well as equipment for food and beverage processing. It is used to decontaminate animal facilities. It is also employed in the health care industries to decontaminate rooms, pass-throughs, isolators and also as a sterilant for product and component sterilization. What's more, as an oxidizing agent, it is extensively used to bleach, deodorize, and detoxify a wide variety of materials, including cellulose, paper-pulp, flour, leather, fats and oils, and textiles. Approximately 4 to 5 million pounds of chlorine dioxide are used daily.

Chemical Properties

Pure chlorine dioxide is an unstable gas and therefore is generated as needed. Although chlorine dioxide has "chlorine" in its name, its chemistry is radically different from that of chlorine. When reacting with other substances, it is weaker and more selective. For example, it does not react with ammonia or most organic compounds. Chlorine dioxide oxygenates products rather than chlorinating them. Therefore, unlike chlorine, chlorine dioxide does not produce environmentally undesirable organic compounds containing chlorine.

Chemical Formula: ClO₂

Molecular Weight: 67.45 g/mole

Melting Point (°C): -59

Boiling Point (°C): +11

Density: 2.4 times that of air

Antimicrobial Properties / Mode of Action

Chlorine dioxide (ClO₂) acts as an oxidizing agent and reacts with several cellular constituents, including the cell membrane of microbes. By "stealing" electrons from them (oxidation), it breaks their molecular bonds, resulting in the death of the organism by the break up of the cell. Since chlorine Dioxide alters the proteins involved in the structure of microorganisms, the enzymatic function is broken, causing very rapid bacterial kills. The potency of chlorine dioxide is attributable to the simultaneous, oxidative attack on many proteins thereby preventing the cells from mutating to a resistant form. Additionally, because of the lower reactivity of chlorine Dioxide, its antimicrobial action is retained longer in the presence of organic matter.

Sporal vs bacterial inactivation?

The difference between sporal and bacterial inactivation can be likened to the difference between sterilization and disinfection. For example, for a chemical agent to be classified as a sterilant, it must be demonstrated to have sporicidal activity. For this reason, spores are commonly used as a challenge in sterilization process development. Disinfection, on the other hand, does not require the complete inactivation of microbial life and in some cases disinfectant claims may be substantiated by demonstrating bacterial inactivation capability.

The following is a simplistic explanation of the differences between a spore and a bacterium. Bacterial endospores are one of the most persistent forms of microbial life and typically require aggressive inactivation procedures. Vegetative bacteria are generally much more easily inactivated than are bacterial endospores. This is primarily because the sensitive areas of bacteria are easily contacted by chemosterilizing agents. The spore, however, has a more complex structure than the vegetative bacterial cell. Its sensitive material is contained within a core and that core is surrounded by a cortex and spore coats. These coats tend to act as a permeability barrier to the entry of chlorine dioxide and other compounds (Knapp).

Environmental Impact

Chlorine dioxide's special properties make it an ideal choice to meet the challenges of today's environmentally concerned world. Actually, chlorine dioxide is an environmentally preferred alternative to elemental chlorine. When chlorine reacts with organic matter, undesirable pollutants such as dioxins and bio-accumulative toxic substances are produced. Thus, the EPA supports the substitution of chlorine dioxide for chlorine because it greatly reduces the production of these pollutants. It is a perfect replacement for chlorine, providing all of chlorine's benefits without any of its weaknesses and detriments. Most importantly, chlorine dioxide does not chlorinate organic material, resulting in significant decreases in trihalomethanes (THMs), haloacetic acids (HAAs) and other chlorinated organic compounds. This is particularly important in the primary use for chlorine Dioxide, which is water disinfection. Other properties of chlorine dioxide make it more effective than chlorine, enabling a lower dose and resulting in a lower environmental impact.

What about liquid? Will CD get into it? Will it kill organisms in it? What is the D-value?

In order to maximize process reproducibility and minimize materials effects when using the ClorDiSys Sterilization Systems, it is best to avoid pools or puddles of liquid water. However, if small amounts of liquid are present the efficacy of chlorine dioxide is not affected. The reason that small amounts of water will not impact sterilization efficacy is that chlorine dioxide is readily soluble in water. The partition coefficient ($C_{ClO_2(H_2O)}/C_{ClO_2(air)}$) of chlorine dioxide at 22 °C and 101 kPa is about 38 (Masschelein). And provided that the quantity of water is small the gas concentration in the water reaches equilibrium quickly.

In any case, the final concentration of chlorine dioxide in the water will be higher than the concentration in the gaseous environment. Furthermore the activity of chlorine dioxide in water is even greater than its activity in the gaseous phase. Its bactericidal, virucidal and sporicidal properties in water have been demonstrated at minimum concentrations of 0.20-0.25 mg/L (aq) with temperature dependent D-values for common water contaminants in the range of 16 to 40 seconds at 30 to 20 °C. For gaseous applications, D-values at 20°C for common indicator organisms are 14-45 seconds at 20-10 mg/L (gas).

Process Advantages:

- Short Cycles
- It is a true gas that distributes rapidly
- Biocidal at low concentration and ambient temperature
- No liquids
- Process effectiveness independent of dew point and condensation
- Extremely low residuals
- Rapid aeration (low-use concentration and minimal adsorption)

CD Antimicrobial Spectrum of Activity:

Vegetative Bacteria:

- Staphylococcus aureus
- Pseudomonas aeruginosa
- Salmonella cholerasuis
- Mycobacterium smegmatis

(* CD Indicator Organism)

Bacterial Spores:

- Bacillus subtilis *
- Bacillus stearothermophilus
- Bacillus pumilus
- Clostridium sporogenes

Fungi:

- Aspergillus niger
- Trychophyton mentagrophytes
- Candida albicans

Viruses:

- Polio Type II (non-lipid)
- Herpes simplex Type I (lipid)
- Parvo Virus

D-Value vs. CD Concentration:

<u>CD Concentration (mg/L)</u>	<u>D-Value (minutes)</u>
10	0.75
20	0.27
30	0.12

REFERENCES:

Jeng, D.K. and Woodworth, A.G. (1990) chlorine Dioxide Gas Sterilization Under Square Wave Conditions. *Applied and Environmental Microbiology*. 56(2):514-519.

Knapp J. E. and Battisti D. L. (2000) *chlorine Dioxide*. In: Block, S. ed. Disinfection, Sterilization, and Preservation. 5th edition. Lippencott, Williams and Wilkens. Philadelphia, PA. Chapter 11., p. 215-227.

Masschelein, W.J. and Rice, R.G. editors. (1979) chlorine Dioxide Chemistry and Environmental Impact of Oxychlorine Compounds. Ann Arbor Science Publishers, Inc., Ann Arbor, MI. p.98, 111-145.

Rosenblatt, D.H., Rosenblatt, A.A. and Knapp, J.E. (1987) Use of chlorine Dioxide Gas as a Chemosterilizing Agent. U.S. Patent 4,681,739.

CSI ClorDiSys Solutions, Inc.

P.O. Box 138, Lebanon, NJ 08833-0138 Tel:(908) 236-4100 Fax:(908) 236-2222

www.clordisys.com