

"The Chlorine Dioxide People"

Providing you with gaseous chlorine dioxide solutions for your decontamination and sterilization needs

## Ceiling Plate (Part #: CSI-ROOM014)

### Description

The Ceiling Plate provides a method to connect ClorDiSys Solutions, Inc. Distribution Plate (CSI-ROOM015) to single room. It is mounted in the room on the ceiling or any available place to provide entry into the room from the distribution plate.

This part is used in conjunction with CSI-Room015 which is mounted in a central location.

The ceiling plate is to be securely screwed and sealed to the ceiling to prevent gas from leaking during the decontamination process.



**Ceiling Plate (Part #: CSI-ROOM014)**  
 Unit Dimensions: 13.0" x 13.0"  
 Unit Weight: 13 lbs

