

CSI ClorDiSys Solutions, Inc.

"The Chlorine Dioxide People"

Providing you with gaseous chlorine dioxide solutions for your decontamination and sterilization needs

Under Door Plate (Part #: CSI-ROOM001)

Description



The Under Door Plate provides a method to connect ClorDiSys Solutions, Inc. gas generators to a room to insert the necessary tubing connections under a door, with NO room/wall/door modifications necessary.

The plate allows for quick valved tubing connections and quick connection of electrical connectors to measure and control the RH for either the Mindox or Cloridox systems.

The user only needs one (1) under the door plate since it moves from room to room. This makes the Mindox or Cloridox system installations simple and easy.



Under Door Plate (Part #: CSI-ROOM001)
Unit Dimensions: 8.5" x 8.5" x 24"
Unit Weight: 5 lbs



CSI ClorDiSys Solutions, Inc.

P.O. Box 549, Lebanon, NJ 08833-0549 Tel: (908) 236-4100 Fax: (908) 236-2222

www.clordisys.com

100-005 Rev B