



# ClorDiSys

*"The Chlorine Dioxide People"*

Providing you with gaseous chlorine dioxide solutions for your decontamination needs

™

## Minidox - L Chlorine Dioxide Gas Generator

### Description

The Minidox - L Decontamination System is a smaller, more economical chlorine dioxide gas generation system designed for use in any chamber under 300 ft<sup>3</sup> (8.5 m<sup>3</sup>) such as an isolator, incubator, HEPA housing, or a Biological Safety Cabinet (BSC). It provides a rapid and highly effective method to decontaminate a target chamber. The Minidox -L Decontamination System is portable in design and easily connected to various targets. It includes a BSC interface plate, carbon-based scrubber and blower, and required tubing.



Minidox-L

Dimensions 23"W x 48"H x 24"D  
(58cm W x 122cm H x 61cm D)  
Weight 200 lbs (91kg)



### Why use the Minidox - L gas generators?

- Replacement for formaldehyde
- Generates a true gas
- Excellent distribution into hard to reach areas
- Does not require tight control of dew point temperatures
- No manual wiping or neutralization
- Does not condense out or breakdown during the process
- Non-carcinogenic
- Non-flammable
- Complete cycle automation
- Portable
- Rapid cycle times
- Quick aeration (can literally be minutes)
- No down time between cycles
- Quick connection to chamber
- Sterilization at ambient temperatures
- Easy to use touch screen operator interface
- Quick consumable replacement

### Utility Requirements at Point of Use

100 - 240 VAC, 15 amp, single phase

 **ClorDiSys**

P.O. Box 549, Lebanon, NJ 08833-0549 Tel:(908) 236-4100 Fax:(908) 236-2222

[www.clordisys.com](http://www.clordisys.com)

08-030 Rev C