



"The Chlorine Dioxide People"



Providing you with gaseous chlorine dioxide solutions for your decontamination and sterilization needs

## Wall Plate

### Description:

The Room Wall Plate provides a method to connect ClorDiSys Solutions, Inc gas generators to a room for quick, easy connections. It is mounted through a studded or block wall to the room to be decontaminated.

The plate allows for quick valved tubing connections and quick connection of electrical connectors to measure and control the RH to either all chlorine dioxide gas generators.



Product Number	Description
PLAT-WALL/B/RH	For block room construction, RH Probe included
PLAT-WALL/B	For block room construction, RH Probe not included
PLAT-WALL/S/RH	For studded room construction, RH Probe included
PLAT-WALL/S	For studded room construction, RH Probe not included



P.O. Box 549, Lebanon, NJ 08833-0549 Tel: (908) 236-4100 Fax: (908) 236-2222

[www.clordisys.com](http://www.clordisys.com)

100-009 Rev E