



"The Chlorine Dioxide People"

Providing you with gaseous chlorine dioxide solutions for your decontamination and sterilization needs

Scrubber Systems

Description:

ClorDiSys Solutions offers a broad range of scrubbers if you have a unique circumstance where direct discharge to the atmosphere is not acceptable. Both carbon and sodium thiosulfate scrubbing technologies are available. We have a few standard units available for applications

such as isolators, Biological Safety Cabinets (BSC's), Sterilizers, Small Pass-through rooms, and pressure relief systems. Larger scrubbers for large rooms and building decontamination are custom designed for the particular application.



Room Scrubber



BSC Scrubber

Product Number	Application	Dimensions
SCRSYS/BSC-001	BSC's, Small Rooms, Isolators, Pass-through Chambers	43" H x 13" W x 18" D (33 cm x 109 cm x 46 cm)
SCRSYS-ROOM/100-001	Room Scrubber	47" H x 22" W x 32" D (56 cm x 119 cm x 81 cm)



P.O. Box 549, Lebanon, NJ 08833-0549 Tel:(908) 236-4100 Fax:(908) 236-2222

www.clordisys.com

100-008 Rev 10-2014